

NEW



VAC

**Vacuum Atmospheres
Company**

In Box Solvent Dryer

HIGHLIGHTS:


- As configured dries up to 75 ml of solvent in minutes.
- Karl Fischer test results indicate <10 PPM in THF.
- Ability to dry just enough solvent needed for a specific reaction.
- Easy snap in and out drying columns.
- User replaceable drying media.
- Convenient in box design eliminates inadvertent exposure to air.
- Does not require a vacuum pump or vacuum lines.
- Does not require an inert gas source.
- Easily mounts in any make/model glove box.



www.vac-atm.com
info@vac-atm.com

Main Offices and Manufacturing Facilities
4652 W. Rosecrans Ave,
Hawthorne, CA 90250 (310) 644-0255
East Coast Office (978) 561-1763

Moisture Removal vs Time



| | |
|------------|---|
| VAC 108442 | System Assembly, In Box Solvent Dryer |
| VAC 108507 | Cartridge Assy INBSP, Tetrahydrofuran (THF), CAS 109-99-9 |
| VAC 108508 | Cartridge Assy INBSP, Hexane, CAS 110-54-3 |
| VAC 108509 | Cartridge Assy INBSP, Dichloromethane, CAS 75-09-2 |
| VAC 108510 | Cartridge Assy INBSP, Pentane, CAS 109-66-0 |
| VAC 108511 | Cartridge Assy INBSP, Methanol, CAS 67-56-1 |
| VAC 108512 | Cartridge Assy INBSP, Diethyl Ether, CAS 60-29-7 |
| VAC 108513 | Cartridge Assy INBSP, Ethylene Glycol Dimethyl Ether, CAS 110-71-4 |
| VAC 108514 | Cartridge Assy INBSP, Toluene, CAS 108-88-3 |
| VAC 108515 | Cartridge Assy INBSP, Triethylamine, CAS 121-44-8 |
| VAC 108516 | Cartridge Assy INBSP, Acetone, CAS 67-64-1 |
| VAC 108517 | Cartridge Assy INBSP, Acetonitrile, CAS 75-05-8 |
| VAC 108518 | Cartridge Assy INBSP, Benzene, CAS 71-43-2 |
| VAC 108519 | Cartridge Assy INBSP, Diisopropyl Ethyl Amine, CAS 7078-68-5 |
| VAC 108520 | Cartridge Assy INBSP, Dimethylformamide (DMF), CAS 68-12-2 |
| VAC 108521 | Cartridge Assy INBSP, Dimethyl-Sulfoxide CAS 67-68-5 |
| VAC 108522 | Cartridge Assy INBSP, Chloroform, CAS 67-66-3 |
| VAC 108523 | Cartridge Assy INBSP, Heptane, CAS 142-82-5 |
| VAC 108524 | Cartridge Assy INBSP, Methyl Tert-Butyl Ether (MTBE), CAS 1634-04-4 |
| VAC 108525 | Cartridge Assy INBSP, Benzotrifluoride, CAS 98-08-8 |
| VAC 108526 | Cartridge Assy INBSP, Dioxane, CAS 123-91-1 |
| VAC 108527 | Cartridge Assy INBSP, Pyridine, CAS 110-86-1 |
| VAC 108528 | Cartridge Assy INBSP, Ethanol, CAS 64-17-5 |
| VAC 108529 | Cartridge Assy INBSP, Dodecane, CAS 112-40-3 |
| VAC 108530 | Cartridge Assy INBSP, Diethylformamide (DEF), CAS 617-84-5 |
| VAC 108531 | Cartridge Assy INBSP, Butanol, CAS 71-36-3 |
| VAC 108532 | Cartridge Assy INBSP, Cyclohexane, CAS 110-82-7 |
| VAC 108533 | Cartridge Assy INBSP, Chlorobenzene, CAS 108-90-7 |
| VAC 108534 | Cartridge Assy INBSP, Propylene Carbonate, CAS 108-32-7 |
| VAC 108535 | Cartridge Assy INBSP, Dimethylacetamide (DMAC), CAS 127-19-5 |
| VAC 108536 | Cartridge Assy INBSP, Hexamethyldisiloxane, CAS 107-46-0 |
| VAC 108537 | Cartridge Assy INBSP, Dichloroethane, CAS 107-06-2 |

SETTING THE STANDARD FOR OVER 50 YEARS

